ENGLISH

Indication:

Glycine is a powder indicated for oral prophylaxis procedures, assisting in the removal of supragingival and subgingival dental plaque, as well as extrinsic stains. It may be used in patients with orthodontic brackets and in individuals in the maintenance phase of periodontal therapy, provided that the gingival pockets do not exceed 5 mm in depth. It must be applied by means of air-polishing equipment specifically intended for this purpose.

Intended user:

Duly qualified dental professionals.

Target population:

Patients undergoing oral prophylaxis procedures, including individuals with fixed orthodontic appliances, patients in periodontal maintenance with gingival pockets of up to 5 mm, and people who require mechanical removal of biofilm and extrinsic stains by means of glycine air polishing.

Fundamentals of its operation and action:

Glycine is a product with fine abrasive particles that, through an air-polishing system, performs oral prophylaxis with removal of dental plaque.

Contraindication:

- Patients with asthma, bronchitis and respiratory problems, as the fine powder mist may aggravate their clinical condition.
- Immunosuppressed patients require great caution in the use of this product, since treating deep periodontal pockets may cause bacteremia.
- Discontinue use of the product in the event of allergic reactions.
- Use is contraindicated on prostheses whose material may be affected by the powder jet (in patients with restorations made of metal, composite resin or ceramic, glycine may be used).

Instructions for use:

- Apply a lip balm for greater patient comfort, preventing dryness and cracking.
- Shake the glycine bottle and, following the manual of the air-polishing equipment manufacturer, fill the reservoir with powder up to the internal step of the reservoir. Before taking it to the oral cavity, discharge the first seconds of the jet into a sink to make sure that water is being released and that the equipment is working normally.
- The suction device must be in place from the beginning to the completion of the procedure.
- Keep a distance of 3 to 5 mm between the equipment tip and the tooth surface, always angled towards the tooth, in order not to injure the gingival tissues. The jet must be kept on each tooth for no more than 5 seconds.
- The angle of the equipment nozzle in relation to the tooth must be between 30° and 60°.

Presentation:

- Package with 60 g or 120 g.
- Particle size: 25 µm.

Cleaning procedure:

Not applicable. This is a chemical product and therefore must not undergo any type of cleaning procedure.

Sterilization procedure:

Not applicable. As it is a chemical product, glycine must not be sterilized or subjected to any decontamination method.

Warnings and Precautions:

- Only trained and duly qualified dental professionals should use this product, as they have the scientific information necessary for its proper use.
- Wear PPE, including gloves, masks, safety glasses and long-sleeved lab coats when handling the product, as is done with all other dental instruments. Also use protective glasses on the patient. The fine powder may cause irritation; always wear protective glasses, including over prescription glasses.
- Never direct the jet onto soft tissues, as this may cause injuries, and air spray into the bloodstream may cause air embolism and emphysematous changes.
- Always use suction throughout the entire procedure.
- The same air-polishing equipment is commonly used with both sodium bicarbonate and glycine; make sure to clean it thoroughly after the procedures, as bicarbonate residues may remain retained and bicarbonate must not be used in subgingival regions.
- Ingesting food or beverages within the following 2 hours may cause extrinsic staining of the teeth.
- Do not use the product after the expiration date or if the packaging has been damaged.
- Shake the bottle before use.
- Keep out of the reach of children.
- Avoid contact with the eyes. In case of contact, rinse immediately with plenty of water.
- LOT NUMBER: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.

Discard:

- The disposal of the material must follow the relevant local legislation.

Shelf life:

The product has a shelf life of 2 years and must be used within this period to ensure its efficacy and safety.

Service life:

Single-use product, with a service life restricted to the moment of clinical application. After being removed from the original packaging and used, any remaining quantity in the device must be discarded immediately.

Conservation and storage:

- Always keep the product in its packaging, tightly closed.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- Avoid contact with moisture.
- Store in a ventilated and dry place, protected from direct sunlight and at room temperature.
- The carrier must be informed of the content. The precautions to be taken during transport are shown on the packaging, thus ensuring the protection of the product from shipment to delivery to the customer. Packages must not be left exposed to the weather.

Other additional information:

- 3R is not responsible if the instructions for use are not followed correctly, generating undesired results in the efficacy and safety of the product/patient.

SENSI PRO - GLICINA
SENSI PRO - GLYCINE
SENSI PRO - GLICINA

REV.01

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you within a maximum of 7 business days from the request.
- Always check the website: www.microdont.com.br, for the most up-to-date version of the instructions for use corresponding to the product.
- In the event of adverse events or any incidents related to the use of the product, immediately notify the manufacturer by email: qualidade@3r.ind.br.

Composition:

Glycine.

Silanized silica.

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Medical device	Indicates that it is a medical device.
	Keep dry	Indicates a health product that needs to be protected from moisture.
	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Manufacturer	Indicates the manufacturer of the health product.
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalogue number	Indicates the manufacturer's reference number so that the health product can be identified.
	Keep away from sunlight	Indicates a health product that needs protection from light sources.
	For professional use only	Indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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