ENGLISH**Indication:**

The So – Soft is an absorbent pad made of non-woven material, indicated for use in dental procedures that require moisture control in the operative field. When placed in the oral cavity, close to the salivary glands, it absorbs saliva and other oral fluids, preventing the build-up of moisture in the working area. Its structure provides high absorption capacity, preserves tissue integrity and minimises the risk of microfibre shedding.

3R Indústria e Comércio Ltda. is not responsible for improper use of the product, especially when used for purposes other than those for which it was designed.

Intended user:

Duly qualified dental professionals.

Target population:

Patients undergoing dental procedures that require moisture control in the operative field. The product is intended for patients of all age groups, subject to the assessment and indication of the qualified professional.

Fundamentals of its operation and action:

The So – Soft works by capillary action, combined with the high porosity of its non-woven structure, which allows saliva and other oral fluids to be absorbed and retained quickly. When placed close to the salivary glands or in areas with excess moisture, it provides passive control of the environment, helping to keep the operative field dry and suitable for the clinical procedure.

Applications:

Moisture control and absorption of saliva and other oral fluids during dental procedures that require a dry operative field.

Contraindications:

- Do not use in patients with a known history of hypersensitivity or allergic reaction to any component of the raw material used in the manufacture of the product.
- Do not use in deep surgical cavities, in areas communicating with the airways, or in situations where there is a risk of retention, aspiration or accidental swallowing of the product.
- Do not use as a substitute for sterile gauze in surgical procedures requiring strict aseptic control.
- Do not reuse the product.

Instructions for Use:

- Remove from the packaging only the required quantity of So – Soft, avoiding contact with the remaining units in order to prevent cross-contamination.
- If necessary, adjust the size of the pad according to the application area.
- It is recommended to shape the pad into a cylindrical form (roll type), which favours adaptation to the salivary regions, such as the areas near the sublingual, submandibular or parotid glands.
- Using clinical tweezers, carefully insert the pad into the oral cavity, positioning it in the desired area.
- If the pad reaches its maximum absorption capacity, remove it with the tweezers and replace it with a new unit.
- Discard the product after use.

Form of presentation:

- Non-sterile product.
- Single use (disposable).

SO – SOFT - COTTON

REV. 00

- Available sizes: 2" x 2" (5,10 cm x 5,10 cm) and 4" x 4" (10,20 cm x 10,20 cm).
- Pack containing 100 or 200 units.

Cleaning Procedure:

- Not applicable. The product is single use and therefore must not be cleaned or reused.

Sterilization Procedure:

- Not applicable. The product is supplied non-sterile and must not be sterilized. Discard after use in accordance with local regulations for dental waste.

Warnings and Precautions:

- The product must be used exclusively by dental surgeons, who have the technical knowledge required for its proper use.
- During the procedure, the professional must wear all personal protective equipment, including protective eyewear.
- Before use, inspect that the packaging is intact and that the So – Soft pad is in perfect condition for use.
- Protective eyewear must also be worn by the patient.
- Non-sterile product.
- Single-use product: do not reuse.
- Keep the product out of the reach of children.
- BATCH NUMBER: it is important to keep the batch number, identified on the product label, so that it can be reported to the manufacturer in the event of any occurrence.

Conservation and storage:

- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.
- The operator must be informed of the content. The precautions to be taken during transport are shown on the packaging, ensuring protection of the product from dispatch to delivery to the customer. Packages must not be exposed to the elements.

Discard:

- Disposal of the material must follow the relevant local legislation.

Shelf life:

- Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Service life:

- Not applicable. The product is for single use only and must be discarded after use.

Additional information:

- If necessary, 3R provides the instructions for use in printed form, at no additional charge, upon request by e-mail (qualidade@3r.ind.br) or by phone (+55 11 5525-0590 / +55 0800-773-5525), and we will respond within a maximum of 7 business days from the request.
- Always check the website www.microdont.com.br for the most up-to-date version of the instructions for use corresponding to the product.

Composition:

SO – SOFT - COTTON

REV. 00

100% cotton.

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Do not reuse	Indicates an intended health product for use or for use on a single patient during a single procedure.
	Keep dry	Indicates a health product that needs to be protected against moisture.
	Medical device	Indicates that it is a medical device.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the batch or batch can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
Rx Only	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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