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ENGLISH

Indication:

Indicated to couple and support dental clinical mirror, providing the dental surgeon with better indirect visualization of intraoral structures during exams and clinical procedures.

Intended user:

Duly qualified dental professionals.

Target population:

Patients who require dental evaluation or treatment with the aid of indirect visualization of intraoral structures, at the discretion of the dental surgeon, without specific restriction as to age or gender.

Fundamentals of its operation and action:

The product acts as a mechanical support for dental clinical mirrors, allowing them to be fixed by means of a standard thread system. Its operation allows the professional to properly position the mirror to obtain an indirect view of intraoral structures, especially in regions with difficult direct visual access, in addition to helping to remove soft tissues and reflect light during clinical examination.

Contraindications:

There are no specific contraindications for the use of the product, as long as it is used according to the instructions for use and by a duly qualified professional.

Instructions for use:

- Before use, make sure that the product is clean, whole, and properly sterilized.
- Attach the clinical mirror to the handle by means of the screw system, ensuring firm and stable fixation.
- Position the assembly (handle + mirror) in the oral cavity, adjusting the angle as necessary to obtain an indirect view of the intraoral structures.
- Use the instrument delicately, which can also help to remove soft tissues and reflect light during clinical examination.
- After use, detach the clinical mirror from the handle to carry out the cleaning, disinfection and sterilization procedures as recommended.

Form of presentation:

- Non-sterile product;
- Package containing 1 (one) unit.

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REV.01

Manual Cleaning Procedure

If manual cleaning is the only option available, proceed as follows:

- Use an area specifically reserved for cleaning;
- Rinse the device under running water and keep it immersed in a bath of neutral pH cleaning solution (detergent enzymatic), following the instructions of the agent manufacturer.
- After soaking and keeping it immersed, brush well away from the body using the neutral cleaning agent (enzymatic detergent). Care should be taken to avoid spreading contaminants by spraying during the brushing process. Do not use wire brushes.
- Special care should be taken to clean crevices and other hard-to-reach areas.
- Visually inspect to confirm debris removal.
- Repeat the cycle if necessary.
- Rinse thoroughly under running water until visibly clean.
- Dry the device using a clean compressed air, or paper towel.

Ultrasonic Cleaning Procedure

For cleaning by means of ultrasound, it is recommended to proceed as follows:

- Prepare a pH-neutral cleaning solution (enzymatic detergent), following the instructions of the agent manufacturer for the correct concentration, exposure time, temperature and water quality.
- Place in the proper location of the ultrasonic equipment, following the equipment manufacturer's instructions for use.
- Completely submerge the device in the cleaning solution and program the equipment for at least fifteen (15) minutes.
- Rinse thoroughly under running water until visibly clean.
- Dry the device using a clean compressed air, or paper towel

Sterilization Requirements

Sterilization in a saturated steam autoclave under pressure is recommended, at 121°C for 20 minutes, or the following parameters should be followed

of sterilization in accordance with the relevant policies/regulations of the country where the device is being used.

Warnings and Precautions:

- Do not use the product if there is damage or failures in the thread system that prevent the proper fixation of the clinical mirror.
- Make sure that the clinical mirror is correctly coupled to the handle, ensuring a firm hold before use.
- The product must be subjected to proper cleaning, disinfection and sterilization processes before each use.
- Inappropriate use of the product may compromise intraoral visualization and the safety of the clinical procedure.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.

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REV.01

Discard:

- The disposal of the material must follow the relevant local legislation.

Shelf life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored as per the conditions

indicated, it maintains its properties indefinitely.

Service life:

The service life of the product may vary according to the conditions of use, frequency of sterilization, storage and maintenance. It is recommended that the product be regularly inspected for signs of wear, deformation or corrosion, and replaced whenever it presents any alteration that may compromise its performance or safety.

Conservation and storage:

- Keep the packaging in a dry place, away from light and away from the action of bad weather.
- Keep products in their original packaging until the moment of use

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you within a maximum of 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

AISI 303 Stainless Steel Mirror Handle

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a medical device that has not undergone a sterilization process.
	Medical device	It indicates that it is a medical device.
	See instructions for use	Indicates the user's need to consult the instructions for use.
	Date of manufacture	Indicates the date the medical device was manufactured.
	Manufacturer	Indicates the manufacturer of the medical.
	Código do lote	Indica el código de lote del fabricante para que se pueda identificar el lote o lote.
	Número de catálogo	Indica el número de referencia del fabricante para que se pueda identificar el producto para la salud.
	Solo para uso profesional	Indica que el dispositivo médico debe ser utilizado exclusivamente por un profesional de la salud calificado.

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REV.01

Rx Only	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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